

fascia releaZer®



MG 850

beurer  
active

## MG 850 Fascia-ReleaZer® and Deep-ReleaZer Fascia massager

german  
engineering



Muscle  
fasciae tissue

Breaks down fascia agglutinations



regeneration

Significantly increases blood flow



attach-  
ments

3

3 attachments with 3 different degrees of hardness

Well-known application from physio practice

### Fascia-ReleaZer®:

Significant increase in circulation

Rejuvenates and regenerates the muscle fasciae tissue

Increases suppleness and mobility significantly

Deeply working vibration massage - 2 programmes

4 treatment surfaces

Battery life up to 11 hours per charge

Charging time approx. 2.5 hours

### Deep-ReleaZer attachments:

- Ideal for back and neck treatment as well as for hip and thigh

- Relieves tension in the back area

- Loosens muscles

- 3 attachments with 3 different degrees of hardness for an individual treatment

- Easy to attach with hook-and-loop fasteners

Includes application DVD

Product measurements (W/H/D): 60 x 5 x 3.7 cm

Product weight: approx. 600 g

Mains voltage:

Input: 100-240V~, 50/60Hz, 160mA

Output: 5V --- 1000mA

3 year guarantee

Sales unit: 2 / Shipping carton: -

EAN no.: 4211125 64038 5

Item no.: 640.38




640.38\_MG850\_2019-05-14\_01\_DS\_EN\_BEU Subject to errors and changes

PhysioLine



**Fascia-ReleaZer®**  
incl. Deep-ReleaZer attachments

 german  
engineering



All muscle fibres are surrounded by fasciae. Regular use of the Fascia-ReleaZer® breaks down fasciae agglutinations and actively aids fitness, performance and personal well-being. Tensions can be relieved easily thanks to the attachments.

**Fascia-ReleaZer®**

- Rejuvenates and regenerates the muscle fasciae tissue
- Significant increase in circulation and cellulose regeneration
- Increases suppleness and mobility significantly

**Deep-ReleaZer attachments:**

- Ideal for back and neck
- Relieves tension in the back area